

### Characteristics



<b>Temperature Range</b>	-10/-60 °C
<b>Capacity</b>	300 lt (10,60 cu ft)
<b>External Dimensions (WxDxH)</b>	cm 130,5x73x86,5 <i>Internal Dimensions (WxDxH) cm 111x45x65 (+ 20x45x40);</i>
<b>Weight</b>	Kg 72
<b>Structure and insulation</b>	White coated steel inside and outside, thickness 100 mm
<b>Feet</b>	n.4 feet
<b>Door</b>	n. 1 lid
<b>Key Lock</b>	
<b>Controller</b>	Digital electronic thermostat
<b>Alarms</b>	Acoustic and Visual for temperature deviations
<b>Refrigeration Type</b>	Static
<b>Defrost</b>	Manual

### Technical data

<b>Details</b>	<ul style="list-style-type: none"> <li>Refrigerant type: Completely eco-friendly mixture of refrigerant gases CFC free</li> <li>Mains Voltage 230V/1/50Hz</li> <li>Instant power consumption: 0,56 kW (4 A)</li> <li>Noise level &lt; 50 dB</li> <li>Ambient Temperature: +10/+32 °C; humidity 60%</li> <li>Connector: Schuko</li> </ul>
<b>In accordance with</b>	2014/30/CE, 2014/35/CE and following modifications
<b>Package</b>	Strong Cardboard on wooden pallet Dim. cm. 140x80x110; Kg.82 <i>Wooden crate on wooden pallet available: dim. cm.155x90x112; Kg.133</i>

Images are purely indicative. The features and the accessories of the products may differ according to model and application.

**Swivel castors (of which 2 with brakes) - H +11 cm**

---

**Internal steel wire basket (max n.3)**

dimensions WxDxH cm 40x25x15

---

**Voltage stabilizer**

---

**Frequency 60 Hz**

---

**Voltage 115 V**

---

**FIOCCHETTI SCIENTIFIC S.R.L.**

Headquarter: Via Panagulis, 48 - 42045 Luzzara (RE) - Italy  
VAT No. IT00934960352

Phone +39 0522 976232 - Fax +39 0522 976028

[www.fiocchetti.com](http://www.fiocchetti.com) - Email: [info@fiocchetti.it](mailto:info@fiocchetti.it) - Pec: [fiocchetti@pec.it](mailto:fiocchetti@pec.it)